

Healthcare Access Tablet 7”

The Healthcare Access Tablet (HAT) is a standalone device with an interactive touchscreen. With built-in Wi-Fi (and optional 4G LTE), it requires only a power connection for set up and automatically connects to the network from virtually anywhere.

- 7” touch screen LCD display
- Built-in Wi-Fi, optional 4G LTE
- Built-in camera for video conferencing
- Wireless medical device connectivity
- Support for wired devices
- Built-in speakers



Healthcare Access Tablet 7

Technical Specification	
Screen Dimensions	5" (H) x 7" (W) x 4" (D)
Weight	Approximately 2 lbs.
Display	7" touchscreen (1024 x 600), HDMI
Camera	1080p max
Network	Wi-Fi, 4G LTE optional
Bluetooth Connectivity	Bluetooth LE (device communication only)
Speaker	2 x 2W
Processor	ARM based Quadcore processor
Power Supply	100-240VAC 50/60Hz 12VDC
Ports	2x USB (reserved for medical devices)
Included in the Box	HAT, Power supply

- Low Cost
- All-in-one tablet with built-in Wi-Fi
- Built-in camera/microphone for instant video-calling
- Simple device software management
- Light weight - weighs about 2 lbs.
- General examination camera
- Built-in stethoscope hardware

*FDA Class II clearance pending

Vital Sign Collection: A wide variety of wired and wireless medical devices can connect to the Tablet to collect the patient's vital signs.

Assessment Questionnaires (Optional): Using our Healthcare Access Server, a series of questions and responses can be created to evaluate the patient's health status.

Video Education (Optional): The patient can be assigned multimedia content that can be viewed directly on the Tablet to improve comprehension and participation in their health care plan.