

Healthcare Access Case

The Healthcare Access Tablet (HAT) is a standalone device with an interactive touchscreen, mounted in a sleek case. With built-in Wi-Fi (and optional 4G LTE), it requires only a power connection for set up and automatically connects to the network from virtually anywhere.

- 15" touch screen LCD display
- Built-in Wi-Fi, optional 4G LTE
- Built-in camera for video conferencing
- Wireless medical device connectivity
- Support for wired devices



Healthcare Access Case

Technical Specification	
Dimensions	17" (H) x 13" (W) x 4" (D)
Weight	Approximately 15.2 lbs.
Display	15" touchscreen (800 x 600), HDMI
Camera	1080p max
Network	Wi-Fi, 4G LTE optional
Bluetooth Connectivity	Bluetooth LE (device communication only)
Speaker	12W
Processor	ARM based Quadcore processor
Power Supply	100-240VAC 50/60Hz 12VDC
Ports	4x USB (reserved for medical devices)
Included in the Box	HAT with Case, power supply, Pulse Oximeter, Blood Pressure meter with cuff, thermometer, stethoscope and general examination camera, battery, keyboard and mouse

- Completely portable solution
- All-in-one tablet with built-in wireless modem
- Built-in camera/microphone for instant video-calling
- Simple device software management
- General examination camera
- Built-in stethoscope hardware

*FDA Class II clearance pending

Vital Sign Collection: A wide variety of wired and wireless medical devices can connect to the Tablet to collect the patient's vital signs.

Assessment Questionnaires (Optional): Using our Healthcare Access Server, a series of questions and responses can be created to evaluate the patient's health status.

Video Education (Optional): The patient can be assigned multimedia content that can be viewed directly on the Tablet to improve comprehension and participation in their health care plan.