

Healthcare Access Wall Mount System

The Healthcare Access Wall Unit is a standalone device with an interactive touchscreen. With built-in Wi-Fi (and optional 4G LTE), it requires only a power connection for set up and automatically connects to the network from virtually anywhere.

- 24" touch HD display
- Built-in Wi-Fi, optional 4G LTE
- Built-in camera for video conferencing
- Wireless medical device connectivity
- Support for wired devices
- Built-in speakers



Healthcare Access Wall Mount System

Technical Specification	
Dimensions	45" (H) x 33" (W) x 9" (D)
Weight	Approximately 62 lbs. without monitor
Display	24" touchscreen (1920 x 1080), HDMI
Camera	1080p max
Network	Wi-Fi, 4G LTE optional
Bluetooth Connectivity	Bluetooth LE (device communication only)
Speaker	12W x2
Processor	ARM based Quadcore processor
Power Supply	100-240VAC 50/60Hz 12VDC
Ports	6x USB (reserved for medical devices)
Included in the Box	Wall Unit, Power supply

- All-in-one tablet with built-in wireless modem
- Built-in camera/microphone for instant video-calling
- Simple device management
- General examination camera
- Built-in stethoscope hardware

*FDA Class II clearance pending

Vital Sign Collection: A wide variety of wired and wireless medical devices can connect to the Tablet to collect the patient's vital signs.

Assessment Questionnaires (Optional): Using our Healthcare Access Server, a series of questions and responses can be created to evaluate the patient's health status.

Video Education (Optional): The patient can be assigned multimedia content that can be viewed directly on the Tablet to improve comprehension and participation in their health care plan.